


GALACTOGRAPHY_____ **NODULE LOCALIZATION**



**26th Year
Manufacturing
Since 1993**





M.D.L. Srl

Via Tavani 1/a 23014 Delebio (So) Italy

Tel +39 0342 682130 Fax +39 0342 691316

email : info@mdlsrl.com web : www.mdlsrl.com



INDEX

MARK	06
AUTO MARK	07
REPO MARK	08
REPO MARK D	09
GALASET	10
GALASET DUO	11
GALASET 6	12
Tables	13

MARK

Non-repositionable needle for pre-operative localization of mammary nodules

Single use needle made of a non-repositionable stainless steel wire for pre-operative localization of non-palpable mammary nodules

Centimeter markings for easy determination of insertion depth

Sliding stopper

Two markings on wire to indicate hook position


Echogenic markers

Lancet tip available in 3 versions

Transparent luer cone with insertion funnel for easy wire loading



Comes with a **clamp** to fix the wire to the skin



REF: PNxxyy

pp.13

REF: PNXxxyy

pp.13



REF: PNZxxyy

pp.13



AUTO MARK

*Semiautomatic non-repositionable
needle for pre-operative localization
of mammary nodules*



REF: PNAxxyy

pp.13

REF: PNAxxyy

pp.14

REF: PNAZxxyy

pp.14




0426

REPO MARK

*Repositionable needle for
pre-operative localization of
mammary nodules*

REF: PRxxyy

pp.14



**Repositionable SINGLE
HOOK needle for
pre-operative localization
of non-palpable mammary
nodules**



Repositionable needle for pre-operative localization of mammary nodules Double Hook

pp.14

Centimeter markings
for easy evaluation of
insertion depth

Clamp to anchor wire to skin

High flexibility wire

Repositionable DOUBLE HOOK needle for pre-operative localization of non palpable mammary nodules

GALASET

Kit for galactography with rounded blunt end-point needle and luer extension

Thin tube to prevent fluid loss during procedure

Rounded blunt tip

Luer extension

Two tip versions available

25 Gauge

REF: PZXCxxyy

pp.14

30 Gauge

REF: PZXCxxyy

pp.14

Straight tip

REF: PZXDxxyy


pp.14

Kit for galactography procedure provided with rounded blunt end-point needle and luer extension device



GALAKIT DUO

*Kit for galactography consisting of
two dilators with diameter 0.25 e 0.30
mm, 1 ml syringe, Luer extension and
rounded blunt tip needle*



Curved tip

REF: PZYCxyxy

pp.14

Curved or linear tip available
with different codes

Straight tip

REF: PZYDxyxy

pp.14

Kit for galactography procedure
consisting of two dilators having
diameter of 0.25 and 0.30 mm, 1
ml syringe, Luer extension and
rounded blunt tip needle



GALASET6

Kit for galactography consisting of
six dilators 1 ml syringe and luer
extention



Thin tube with luer
extension



Rounded blunt needle tip



Six dilators with handles for rotation



Curved tip

REF: PZZCxyyy

pp.14

Straight tip

REF: PZZDxyyy

pp.14

Kit for galactography
procedure consisting of six
dilators, 1 ml syringe and
luer extension device



ORDER INFORMATION



Mark - pp.06


CODE	Gauge	Length
PN1805	18G	50 mm
PN1807	18G	70 mm
PN1810	18G	100 mm
PN1812	18G	120 mm
PN1815	18G	150 mm
PN2005	20G	50 mm
PN2007	20G	70 mm
PN2010	20G	100 mm
PN2012	20G	120 mm
PN2015	20G	150 mm
PN2105	21G	50 mm
PN2107	21G	70 mm
PN2110	21G	100 mm
PN2112	21G	120 mm
PN2115	21G	150 mm

Mark Z - pp.06

CODE	Gauge	Length
PNZ1805	18G	50 mm
PNZ1807	18G	70 mm
PNZ1810	18G	100 mm
PNZ1812	18G	120 mm
PNZ1815	18G	150 mm
PNZ2005	20G	50 mm
PNZ2007	20G	70 mm
PNZ2010	20G	100 mm
PNZ2012	20G	120 mm
PNZ2015	20G	150 mm

Mark X - pp.06

CODE	Gauge	Length
PNX1805	18G	50 mm
PNX1807	18G	70 mm
PNX1810	18G	100 mm
PNX1812	18G	120 mm
PNX1815	18G	150 mm
PNX2005	20G	50 mm
PNX2007	20G	70 mm
PNX2010	20G	100 mm
PNX2012	20G	120 mm
PNX2015	20G	150 mm



AutoMark - pp.07

CODE	Gauge	Length
PNA1805	18G	50 mm
PNA1807	18G	70 mm
PNA1810	18G	100 mm
PNA1812	18G	120 mm
PNA1815	18G	150 mm
PNA2005	20G	50 mm
PNA2007	20G	70 mm
PNA2010	20G	100 mm
PNA2012	20G	120 mm
PNA2015	20G	150 mm
PNA2105	21G	50 mm
PNA2107	21G	70 mm
PNA2110	21G	100 mm
PNA2112	21G	120 mm
PNA2115	21G	150 mm



AutoMark X - pp.07

CODE	Gauge	Length
PNAX1805	18G	50 mm
PNAX1807	18G	70 mm
PNAX1810	18G	100 mm
PNAX1812	18G	120 mm
PNAX1815	18G	150 mm
PNAX2005	20G	50 mm
PNAX2007	20G	70 mm
PNAX2010	20G	100 mm
PNAX2012	20G	120 mm
PNAX2015	20G	150 mm
PNAX2105	21G	50 mm
PNAX2107	21G	70 mm
PNAX2110	21G	100 mm
PNAX2112	21G	120 mm
PNAX2115	21G	150 mm

AutoMark Z - pp.07

CODE	Gauge	Length
PNAZ1805	18G	50 mm
PNAZ1807	18G	70 mm
PNAZ1810	18G	100 mm
PNAZ1812	18G	120 mm
PNAZ1815	18G	150 mm
PNAZ2005	20G	50 mm
PNAZ2007	20G	70 mm
PNAZ2010	20G	100 mm
PNAZ2012	20G	120 mm
PNAZ2015	20G	150 mm
PNAZ2105	21G	50 mm
PNAZ2107	21G	70 mm
PNAZ2110	21G	100 mm
PNAZ2112	21G	120 mm
PNAZ2115	21G	150 mm



RepoMark - pp.08

CODE	Gauge	Length
PR2005	20G	50 mm
PR2007	20G	70 mm
PR2009	20G	90 mm
PR2010	20G	100 mm
PR2012	20G	120 mm
PR2015	20G	150 mm

RepoMark D - pp.09

CODE	Gauge	Length
PRD2005	20G	50 mm
PRD2007	20G	70 mm
PRD2009	20G	90 mm
PRD2010	20G	100 mm
PRD2012	20G	120 mm
PRD2015	20G	150 mm

Galaset - pp.10

CODE	Gauge	Length	Tip
PZXD2515	25G	150 mm	STRAIGHT
PZXC2515	25G	150 mm	CURVED
PZXD3015	30G	150 mm	STRAIGHT
PZXC3015	30G	150 mm	CURVED

GalaKit Duo- pp.11

CODE	Gauge	Length	Tip
PZYC2515	25G	150 mm	CURVED
PZYD3015	30G	150 mm	STRAIGHT
PZYC3015	30G	150 mm	CURVED

GalaKit 6 - pp.12

CODE	Gauge	Length	Tip
PZZD3015	30G	150 mm	STRAIGHT
PZZC3015	30G	150 mm	CURVED



M.D.L. Srl
Via Tavani 1/a 23014 Delebio (So) Italy
Tel +39 0342 682130 Fax +39 0342 691316
email : info@mdlsrl.com web : www.mdlsrl.com

